



## Guargum Powder - 5000 CPS/200 Mesh

### Technical Specifications:

<b>Viscosity:</b> 1% solution in cold water, cps (As is basis) 1 Hour 24 Hours Brookfield RVF Viscometer, Spindle No. 4, 20 RPM 20 Deg C. PH Granulometry (%) ON the US mesh 100(150 micron) Thru US mesh 200 (75 micron)
<b>Composition:</b> Moisture % Residue insoluble in Acid % Gum % (by difference)
<b>Heavy Metals:</b> Total Heavy Metals (As Pb)
<b>Typical Microbial:</b> Aerobic Plate Count/g Mold and Yeast Coliform /g (as MPN) E. Coli/1.25 g Salmonella

**Application:** Food manufacturing

**Source:** India

**Packaging:** 25 Kg bags  
1000 kg bags  
1000 lb bags

Custom Packaging Available