



## Guargum Powder - 6000 – 7000 CPS/200 Mesh

### Technical Specifications:

<b>Viscosity:</b>	
1% solution in cold water, cps (As is basis)	
10 minutes	4500 - 5000
2 Hours	6000 - 6500
Brookfield RVF Viscometer, Spindle No. 4, 20 RPM	
20 Deg C.	
PH	
Granulometry (%)	10.0 Max
ON the US mesh 100(150 micron)	90.0 Min
Thru US mesh 200 (75 micron)	
<b>Composition:</b>	
Moisture %	10.0 Max
Residue insoluble in Acid %	3.0 Max
Gum % (by difference)	87.0 Min
<b>Heavy Metals:</b>	
Total Heavy Metals (As Pb)	<20 ppm
<b>Typical Microbial:</b>	
Aerobic Plate Count/g	
Mold and Yeast	
Coliform /g (as MPN)	
E. Coli/1.25 g	
Salmonella	

**Application:** Oil & Gas – Drilling

**Source:** India

**Packaging:** 25 Kg bags  
1000 kg bags  
1000 lb bags

Custom Packaging Available